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Indications

Indicated for primary implantation for the visual correction of aphakia in adult patients in whom the cataractous lens has been removed. The lens is intended for placement in the capsular bag.

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Contraindications

There are no known contraindications for the enVista IOLs.

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Warnings

Physicians considering lens implantation under any of the following circumstances should weigh the potential risk/benefit ratio:

1. Recurrent severe anterior or posterior segment inflammation or uveitis.
2. Patients in whom the intraocular lens may affect the ability to observe, diagnose, or treat posterior segment diseases.
3. Surgical difficulties at the time of cataract extraction, which might increase the potential for complications (e.g., persistent bleeding, significant iris damage, uncontrolled positive pressure, or significant vitreous prolapse or loss).
4. A distorted eye due to previous trauma or developmental defect in which appropriate support of the IOL is not possible.
5. Circumstances that would result in damage

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to the endothelium during implantation.

6. Suspected microbial infection.
7. Patients in whom neither the posterior capsule nor zonules are intact enough to provide support.

Precautions

1. Do not attempt to resterilize the lens as this can produce undesirable side effects.
2. Do not use if product sterility or quality is thought to be compromised due to damaged packaging or signs of leakage (such as the loss of saline storage solution, or the presence of salt crystallization).
3. Do not soak or rinse the intraocular lens with any solution other than sterile balanced salt solution or sterile normal saline.
4. Do not store the lens at a temperature greater than 43°C. DO NOT FREEZE. Do not autoclave the intraocular lens.
5. Do not reuse the lens. It is intended for permanent implantation. If explanted, sterility and proper function cannot be assured.
6. The safety and effectiveness of the enVista IOL have not been substantiated in patients with preexisting ocular conditions and intraoperative complications (see below). Careful preoperative evaluation and sound clinical judgment should be used by the surgeon to decide the benefit/risk ratio before implanting a lens in a patient with one or more of these conditions. Physicians considering lens implantation in such patients should explore the use of alternative methods of aphakic correction and consider lens implantation only if alternatives are deemed unsatisfactory in meeting the needs of the patient.
7. If this product is reprocessed and/or re-used,

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Bausch & Lomb cannot guarantee the functionality, material structure, or cleanliness or sterility of the product. Re-use could lead to illness, infection and/or injury to the patient or user and, in extreme incidents, death. This product is labeled as 'single use' which is defined as a device intended to be used once only for a single patient.

Before Surgery

• Retinal conditions or predisposition to retinal conditions, previous history of, or a predisposition to, retinal detachment or proliferative diabetic retinopathy, in which future treatment may be compromised by implanting this lens.

• Amblyopia

• Clinically severe corneal dystrophy (e.g., Fuchs')

• Rubella, congenital, traumatic or complicated cataracts

• Extremely shallow anterior chamber, not due to swollen cataract

• Recurrent anterior or posterior segment inflammation of unknown etiology, or any disease producing an inflammatory reaction in the eye (e.g. iritis or uveitis).

• Aniridia

• Iris neovascularization

• Glaucoma (uncontrolled or controlled with medication)

• Microphthalmos or macrophthalmos

• Optic nerve atrophy

• Previous corneal transplant

• Pre-existing ocular conditions which may negatively impact stability of the implant.

During Surgery

• Mechanical or surgical manipulation required to enlarge the pupil

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• Uncontrollable positive intraocular pressure

• Complications in which the IOL stability

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could be compromised

8. Patients with preoperative problems such as corneal endothelial disease, abnormal cornea, macular degeneration, retinal degeneration, glaucoma, and chronic drug miosis may not achieve the visual acuity of patients without such problems. The physician must determine the benefits to be derived from lens implantation when such conditions exist.

9. A high level of surgical skill is required for intraocular lens implantation. The surgeon should have observed and/or assisted in numerous implantations and successfully completed one or more courses on intraocular lens implantation before attempting to implant intraocular lenses.

10. As with any surgical procedure, there is risk involved. Potential complications accompanying cataract or implant surgery may include, but are not limited to the following: corneal endothelial damage, infection (endophthalmitis), retinal detachment, vitritis, cystoid macular edema, corneal edema, pupillary block, cyclitic membrane, iris prolapse, hypopyon, transient or persistent glaucoma, and secondary surgical intervention. Secondary surgical interventions include, but are not limited to: lens repositioning, lens replacement, vitreous aspiration or iridectomy for pupillary block, wound leak repair, and retinal detachment repair.

11. Care should be taken to remove all viscoelastic from the anterior and posterior surfaces of the lens.

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Potential complications / Side effects

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enVista IOL's are intended to be folded and placed in the posterior chamber of the eye after the cataractous crystalline lens has been removed.

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The recommended A-constant listed on the lens carton is intended for use with axial length measurements obtained by optical biometry. Use of axial length measurements by other techniques (e.g. Applanation A-scan) will

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- Hoffer KJ. The Hoffer Q formula: a comparison of theoretic and regression formulas, Journal of Cataract and Refractive Surgery Vol. 19, pp. 700-712, 1993; ERRATA, Vol. 20, pp. 677, 1994.
- Holladay JT, Musgrove KH, Prager TC, Lewis JW, Chandler TY, Ruiz RS. A three-part system for refining intraocular lens power calculations. Journal of Cataract and Refractive Surgery, Vol. 14, pp. 17-24, 1988.
- Norrby NES. Unfortunate Discrepancies, Letter to the Editor and Reply by Holladay JT. Journal of Cataract and Refractive Surgery, Vol. 24, pp. 433-434, 1998.
- Olsen T, Olesen H, Thim K, and Corydon L. Prediction of pseudophakic anterior chamber depth with the newer IOL calculation formulas. Journal of Cataract and Refractive Surgery, Vol. 18, pp. 280-285, 1992.
- Retzlaff JA, Sanders DR, Kraff MC. Development of the SRK/T intraocular lens implant power calculation formula. Journal of Cataract and Refractive Surgery, Vol. 16, pp. 333-340, 1990; ERRATA, Vol. 16, pp. 528, 1990.
- Haigis W: The Haigis Formula. In: Intraocular lens power calculations. H. John Shamma (eds), Slack Incorporated, Thorofare, NJ, USA, pp. 39-57, 2004.

normally require a different lens constant. This number is a guideline only and is based on an evaluation of clinical data obtained using the IOL Master.

The physician should determine preoperatively the power of the lens to be implanted. Lens power calculation methods are described in the following references:

- Hoffer KJ. The Hoffer Q formula: a comparison of theoretic and regression formulas, Journal of Cataract and Refractive Surgery Vol. 19, pp. 700-712, 1993; ERRATA, Vol. 20, pp. 677, 1994.
- Holladay JT, Musgrove KH, Prager TC, Lewis JW, Chandler TY, Ruiz RS. A three-part system for refining intraocular lens power calculations. Journal of Cataract and Refractive Surgery, Vol. 14, pp. 17-24, 1988.
- Norrby NES. Unfortunate Discrepancies, Letter to the Editor and Reply by Holladay JT. Journal of Cataract and Refractive Surgery, Vol. 24, pp. 433-434, 1998.
- Olsen T, Olesen H, Thim K, and Corydon L. Prediction of pseudophakic anterior chamber depth with the newer IOL calculation formulas. Journal of Cataract and Refractive Surgery, Vol. 18, pp. 280-285, 1992.
- Retzlaff JA, Sanders DR, Kraff MC. Development of the SRK/T intraocular lens implant power calculation formula. Journal of Cataract and Refractive Surgery, Vol. 16, pp. 333-340, 1990; ERRATA, Vol. 16, pp. 528, 1990.
- Haigis W: The Haigis Formula. In: Intraocular lens power calculations. H. John Shamma (eds), Slack Incorporated, Thorofare, NJ, USA, pp. 39-57, 2004.

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1. Prior to implanting, examine the lens package for type, power, and proper configuration.
2. Open the peel pouch and remove the vial in

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a sterile environment.

3. Remove the lid from the vial.

4. With a pair of smooth forceps, remove the lens from the vial by gently grasping the lens haptic.

5. Rinse the entire lens with sterile balanced salt solution or sterile normal saline.

6. Examine the lens thoroughly to ensure particles have not become attached to it, and examine the lens optical surfaces for other defects.

7. The lens may be soaked in sterile balanced salt solution until ready for implantation.

8. Amvisc, Amvisc Plus, or OcuCoat viscoelastic should be used for lubrication of the delivery system / injector when inserting the lens.

9. Bausch + Lomb recommends using a Bausch + Lomb approved delivery system / injector.

10. There are various surgical procedures that can be utilized, and the surgeon should select a procedure that is appropriate for the patient. Surgeons should verify that appropriate instrumentation is available prior to surgery.

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Description of accessories, medical devices or products that are not medical, but intended for use in combination with the declared medical device

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The following compatible medical devices are used for implantation of the enVista lens:

1. Injector for introducing intraocular lenses, INJ100, Registration certificate № RZN 2020/10023;

2. Injection system Bausch + Lomb BLIS, Registration certificate № RZN 2017/6619;

3. Viscoprotector Amvisc, Amvisc Plus with an insertion device, Registration

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Clinical studies have been conducted on the enVista single-piece IOL (model MX60) and the parent hact X-60 three-piece IOL (model X-60). The results of clinical investigations provide reasonable assurance that the Model MX60 IOL is safe and effective for the visual correction of aphakia following cataract extraction.

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For the enVista MX60 clinical study, all subjects in the safety analysis set were evaluated for IOL glistenings at Form 3 and Form 4 visits. IOL glistenings were evaluated via retroillumination slit lamp examination utilizing a photographic grading scale provided in the protocol. The grading scale consisted of (in order of severity), "none, grade 0 (trace), grade 1, 2, 3, or 4." No glistenings of any grade were reported for any subject at any visit in the clinical study.

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The IOL is packaged in a consumer packaging containing:

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- Sterile vial with saline solution 0.9% and a lens inside in a pouch for transportation – 1 pc.;
- Instructions for use – 1 pc.;
- Implant registration patient card – 1 pc.;
- Self-adhesive labels – 8 pcs.

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Intraocular lens is an implantable medical

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Environmental requirements, procedure and conditions for disposal or destruction

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Lenses must be disposed of in accordance with the hygienic requirements of SanPin 2.1.3684-21 "Sanitary and epidemiological requirements for the maintenance of territories of urban and rural settlements, water bodies, drinking water and drinking water supply, atmospheric air, soil, residential premises, operation of industrial and public premises, organization and implementation of sanitary and anti-epidemic (preventive) measures" adopted in the Russian Federation.

Products that have come into contact with blood and / or other biological fluids should be disposed of and / or destroyed in accordance with SanPin 2.1.3684-21 class B.

Unused products (without contact with blood and / or body fluids), incl. Expired in life, subject to disposal and / or destruction in accordance with SanPin 2.1.3684-21 class A.

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Email: Russia.Claims@bausch.com

In case of any questions related to the quality of the medical device, please contact the Authorized Representative of the Manufacturer in the Russian Federation at the following address:

«Bausch Health» LLC, Pavlovskaya St. 7, bldg. 1, premises 1N, int. ter. Danilovsky Municipal District, Moscow, 115093, Russia, tel./fax: +7 495 510 2879

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