Ophthalmic kit for cataract surgery LUXSMART™ YELLOW intraocular lens	Instructions for use for the Russian Federation version 2
Cutting Edge SAS, France	770, Rue Alfred Nobel, Immeuble le Nobel, 34000 Montpellier, France

BAUSCH + LOMB

Instructions for use

Ophthalmic kit for cataract surgery in the following versions LUXSMART™

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ASPHERIC LENS

Model LuxSmart™ (LUXSMART™ yellow preloaded hydrophobic lens)

LUXSMART™, consisting of: 1. Intraocular lens LUXSMART™ YELLOW 2. Delivery system for intraocular lens implantation (Accuject Pro disposable) 3. Instructions for use

DESCRIPTION: The LuxSmart™ intraocular lens (IOL) is a single-piece ultraviolet (UV) absorbing posterior chamber IOL developed to replace the natural crystalline lens in pseudophakic adult patients. The IOL has aspheric optics. The lens is provided sterile and is designed to be folded prior to insertion in the eye and implantation in the capsular bag. The preloaded intraocular LuxSmart™ lens is supplied with an Accuject Pro Delivery System to provide a touch-free method of delivery of the lens into the capsular bag.

PHYSICAL CHARACTERISTICS:

Lens / Haptic Material	Foldable hydrophobic acrylic copolymer containing a UV absorber and a blue light filtering chromophore.												
Optic type	Aspheric												
Index of Refraction when in the eye at 35°C	1,540												
Spectral transmittance	10 % at 396,7 nanometers for +20,0 Diopter												
A-Constant:													
Optical biometry, SRK/T	119,2												
Applanation A-scan, SRK/T	118,8												
Body Diameter	6,0 mm												
Overall Diameter	11,0 mm												
Haptic Angle	0°												
Powers	+0,00 D to +10,00 D with 1,0 D increment +10,50 D to +34,00 D with 0,5 D increment												
	<table border="1"> <thead> <tr> <th colspan="2">Specification</th> </tr> <tr> <th>Power Range</th> <th>Tolerances</th> </tr> </thead> <tbody> <tr> <td>+00.00 D to +15.00 D</td> <td>± 0.30D</td> </tr> <tr> <td>+15.10 D to +25.00 D</td> <td>± 0.40D</td> </tr> <tr> <td>+25.10 D to +30.00 D</td> <td>± 0.50D</td> </tr> <tr> <td>> +30.00D</td> <td>± 1.00D</td> </tr> </tbody> </table>	Specification		Power Range	Tolerances	+00.00 D to +15.00 D	± 0.30D	+15.10 D to +25.00 D	± 0.40D	+25.10 D to +30.00 D	± 0.50D	> +30.00D	± 1.00D
Specification													
Power Range	Tolerances												
+00.00 D to +15.00 D	± 0.30D												
+15.10 D to +25.00 D	± 0.40D												
+25.10 D to +30.00 D	± 0.50D												
> +30.00D	± 1.00D												

Ophthalmic kit for cataract surgery LUXSMART™ YELLOW intraocular lens	Instructions for use for the Russian Federation version 2
Cutting Edge SAS, France	770, Rue Alfred Nobel, Immeuble le Nobel, 34000 Montpellier, France

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MODE OF ACTION:

LuxSmart™ IOL positioned in the capsular bag acts as a refractive medium to correct aphakia.

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Possible adverse events

As with any surgical procedure, there is risk involved. Potential complications of cataract surgery may include, but are not limited to the following: inflammation/infection potentially due to adverse reaction to device materials; complications usually associated with intraocular lens implantation including: Posterior capsule opacification (secondary cataract), inflammation (iritocyclitis, pupillary/cyclitic membrane, vitreous inflammation, CME, TASS ...), corneal endothelial damage, endophthalmitis, retinal detachment, corneal edema, pupillary block, iris prolapse, hypopyon, glaucoma, wound leak, hyphema, iris atrophy, vitreous prolapse, IOL decentration or tilt, IOL optic deposits or opacities, glistenings and secondary surgical intervention (such as but not limited to wound leak repair, vitreous aspiration, iridectomy, IOL repositioning or replacement and retinal detachment repair, etc.).

Ophthalmic kit for cataract surgery LUXSMART™ YELLOW intraocular lens	Instructions for use for the Russian Federation version 2
Cutting Edge SAS, France	770, Rue Alfred Nobel, Immeuble le Nobel, 34000 Montpellier, France

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Do not implant the lens in the anterior chamber or in the ciliary sulcus.

PRECAUTIONS:

- Store at room temperature. Avoid temperatures below 5°C and above 45°C. Do not freeze.
- The LuxSmart™ IOL and the viscoelastic solution used should be stored for at least 30 min at minimum 20°C before surgery.
- A high level of surgical skill is required to implant intraocular lenses. Surgeon should be fully trained before attempting to implant intraocular lenses.
- The IOL model, its power and expiration date should be verified before opening the protective packaging and before opening the individual sterile pouch/blister.
- Sterility of the device is guaranteed only if the individual sterile pouch/blister has not been opened or damaged.
- Handle the lens carefully to avoid surface or haptics damages.
- The IOL must be carefully rinsed prior to implantation. Do not soak or rinse the IOL with any solution other than sterile balanced salt solution.
- Care should be taken to remove all viscoelastic from the anterior and posterior surfaces of the lens after implantation in the capsular bag.
- Care should be taken to verify that the lens is fully implanted into the capsular bag.
- If a YAG laser posterior capsulotomy is performed, assure that the laser beam is focused slightly behind the posterior capsule.

Ophthalmic kit for cataract surgery LUXSMART™ YELLOW intraocular lens	Instructions for use for the Russian Federation version 2
Cutting Edge SAS, France	770, Rue Alfred Nobel, Immeuble le Nobel, 34000 Montpellier, France

PACKAGING / STERILIZATION:

The LuxSmart™ IOL is individually package preloaded in an Accuject Pro delivery system in a blister. All models are supplied dry and should be opened under sterile conditions. A patient card and self-adhesive labels are supplied in the carton to provide traceability of the IOL. All models are sterilized by Ethylene Oxide.

DIRECTIONS FOR USE:

1. The physician should determine preoperatively the power of the lens to be implanted, using calculation methods. The recommended A-constants listed on the lens carton are guidelines only.
2. Prior to implanting, examine the lens package for model, power, expiration date and proper configuration.
3. Open the carton and inspect the sterilized blister containing the lens and the injection system. Make sure it is not damaged. Verify that the information on the blister label is consistent with the information on the outer box and the self-adhesive labels. Peel open the blister and transfer the device onto the sterile field.
4. Inspect the device and make sure it is intact. The lens should be rinsed thoroughly by flushing balanced salt solution through the tip of the cartridge (**FIG. 1**) and from the end of the loading chamber (**FIG. 2**). Wait for 30 seconds.

Please refer to the enclosed figures to fully prepare the device before implantation, taking care not to touch the IOL.

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illustrated in **FIG. 5**, when the haptic insertions are top right (A) and bottom left (B), you are looking at the anterior side of the lens.


FIG. 1


FIG. 2


FIG. 3


FIG. 4


FIG. 5


FIG. 6

Only viscoelastic solution of low to medium viscosity should be used. Only the injection system provided with the preloaded LuxSmart™ lens should be used.

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Transportation conditions:

Temperature range from -10°C to + 40°C and relative humidity from 85% to 90%.

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Manufacturer

A-Constant (est)

Ophthalmic kit for cataract surgery
LUXSMART™ YELLOW intraocular lens

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Tel.: +7-495-510-28-79.

Specifications of Accuject Pro

Inner diameter at the head of the Cartridge tip	1.38 ± 0,05 mm
Outer tip diameter at the head of the Cartridge tip	1.66 ± 0.05 mm
Angle at the head of the tip	45 ± 2°
Total length	151 ± 0.1 mm
Length of the Loading Chamber	18.5 ± 1.0 mm
Width of the Loading chamber (1)	13.53 ± 0.04 mm
Length of the Handpiece	95.77 ± 0.05 mm

Ophthalmic kit for cataract surgery LUXSMART™ YELLOW intraocular lens	Instructions for use for the Russian Federation version 2
Cutting Edge SAS, France	770, Rue Alfred Nobel, Immeuble le Nobel, 34000 Montpellier, France

Width of the Handpiece (1)	35.025 ± 0.03 mm
Width of the Handpiece (2)	8.046 ± 0.03 mm
Length of the Plunger	106.88 ± 0.05 mm
Width of the Plunger (1)	14.927 ± 0.07 mm
Width of the Plunger (2)	5.483 ± 0.03 mm
Width of the Plunger (3)	1.755 ± 0.03 mm
Length of the Silicone tip	10.05 ± 0.05 mm
Width of the Silicone tip (1)	2.348 ± 0.04 mm
Width of the Silicone tip (2)	3.532 ± 0.04 mm
Length of the Cartridge tip	29.11 ± 0.08 mm
Width of the Cartridge tip (1)	10.231 ± 0.04 mm
Width of the Cartridge tip (2)	5.537 ± 0.04 mm
Pressing force on the plunger Injection force	20N ± 2,2 N

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Ophthalmic kit for cataract surgery LUXGOOD™ YELLOW intraocular lens	Instructions for use for the Russian Federation version 2
Cutting Edge SAS, France	770, Rue Alfred Nobel, Immeuble le Nobel, 34000 Montpellier, France

BAUSCH + LOMB

Instructions for use

Ophthalmic kit for cataract surgery in the following versions LUXGOOD™.

CUTTING EDGE

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770 rue Alfred Nobel, 34000 MONTPELLIER

Tél : (33) (0)5 62 24 64 64 - Fax (33) (0)5 62 24 64 66

SAS au capital de 286 309 Euros - RCS Montpellier 800 783 284

N° TVA/VAT : FR 43 800 783 284

Approved by:

28 JUN. 2022

Sandrine BONACHE
CUTTING EDGE SAS
Quality and Regulatory Affairs Manager

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Ophthalmic kit for cataract surgery LUXGOOD™ YELLOW Intraocular lens	Instructions for use for the Russian Federation version 2
Cutting Edge SAS, France	770, Rue Alfred Nobel, Immeuble le Nobel, 34000 Montpellier, France

ASPHERIC LENS

Model LuxGood™ (LUXGOOD™ yellow preloaded hydrophobic lens)

LUXGOOD™, consisting of: 1. Intraocular lens LUXGOOD™ YELLOW 2. Delivery system for intraocular lens implantation (Accuject Pro disposable) 3. Instructions for use

DESCRIPTION: The LuxGood™ intraocular lens (IOL) is a single-piece ultraviolet (UV) absorbing posterior chamber IOL developed to replace the natural crystalline lens in pseudophakic adult patients. The IOL has aspheric optics. The lens is provided sterile and is designed to be folded prior to insertion in the eye and implantation in the capsular bag. The preloaded intraocular LuxGood™ lens is supplied with an Accuject Pro Delivery System to provide a touch-free method of delivery of the lens into the capsular bag.

PHYSICAL CHARACTERISTICS:

Lens / Haptic Material	Foldable hydrophobic acrylic copolymer containing a UV absorber and a blue light filtering chromophore.												
Optic type	Aspheric												
Index of Refraction when in the eye at 35°C	1,540												
Spectral transmittance	10 % at 396,7 nanometers for +20,0 Diopter												
A-Constant:													
Optical biometry, SRK/T	119,2												
Applanation A-scan, SRK/T	118,8												
Body Diameter	6,0 mm												
Overall Diameter	11,0 mm												
Haptic Angle	0°												
Optical Power	+0,00 D to +10,00 D with 1,0 D increment +10,50 D to +34,00 D with 0,5 D increment												
	<table border="1"> <thead> <tr> <th colspan="2">Specification</th> </tr> <tr> <th>Power Range</th> <th>Tolerances</th> </tr> </thead> <tbody> <tr> <td>+00.00 D to +15.00 D</td> <td>± 0.30D</td> </tr> <tr> <td>+15.10 D to +25.00 D</td> <td>± 0.40D</td> </tr> <tr> <td>+25.10 D to +30.00 D</td> <td>± 0.50D</td> </tr> <tr> <td>> +30.00D</td> <td>± 1.00D</td> </tr> </tbody> </table>	Specification		Power Range	Tolerances	+00.00 D to +15.00 D	± 0.30D	+15.10 D to +25.00 D	± 0.40D	+25.10 D to +30.00 D	± 0.50D	> +30.00D	± 1.00D
Specification													
Power Range	Tolerances												
+00.00 D to +15.00 D	± 0.30D												
+15.10 D to +25.00 D	± 0.40D												
+25.10 D to +30.00 D	± 0.50D												
> +30.00D	± 1.00D												

Ophthalmic kit for cataract surgery LUXGOOD™ YELLOW intraocular lens	Instructions for use for the Russian Federation version 2
Cutting Edge SAS, France	770, Rue Alfred Nobel, Immeuble le Nobel, 34000 Montpellier, France

INDICATIONS:

The posterior chamber lens is indicated for the visual correction of aphakia with or without pre-existing corneal astigmatism secondary to removal of a cataractous lens in adult patients. These devices are intended to be placed in the capsular bag.

MODE OF ACTION:

LuxGood™ IOL positioned in the capsular bag acts as a refractive medium to correct aphakia.

CONTRA-INDICATIONS AND WARNINGS:

Contra-Indications

Physicians considering lens implantation in any of the following circumstances (but not limited to) must weigh the potential risk/benefit ratio: uncontrolled positive pressure (glaucoma), recurrent severe anterior or posterior segment inflammation or uveitis; patients in whom the IOL may affect the ability to observe, diagnose, or treat posterior segment diseases; surgical difficulties increasing the potential for complications (such as but not limited to: persistent bleeding, significant iris damage, uncontrolled positive pressure, or significant vitreous prolapse or loss); distorted eye due to previous trauma or developmental defect in which appropriate support of the IOL is not possible; circumstances that would result in damage to the endothelium during implantation; patients in whom neither the posterior capsule nor zonule are intact enough to provide support; suspected microbial infection.

Possible adverse events

As with any surgical procedure, there is risk involved. Potential complications of cataract surgery may include, but are not limited to the following: inflammation/infection potentially due to adverse reaction to device materials; complications usually associated with intraocular lens implantation including: Posterior capsule opacification (secondary cataract), inflammation (iritis, iridocyclitis, pupillary/cyclitic membrane, vitreous inflammation, CME, TASS ...), corneal endothelial damage, endophthalmitis, retinal detachment, corneal edema, pupillary block, iris prolapse, hypopyon, glaucoma, wound leak, hyphema, iris atrophy, vitreous prolapse, IOL decentration or tilt, IOL optic deposits or opacities, glistenings and secondary surgical intervention (such as but not limited to wound leak repair, vitreous aspiration, iridectomy, IOL repositioning or replacement and retinal detachment repair, etc.).

Ophthalmic kit for cataract surgery LUXGOOD™ YELLOW intraocular lens	Instructions for use for the Russian Federation version 2
Cutting Edge SAS, France	770, Rue Alfred Nobel, Immeuble le Nobel, 34000 Montpellier, France

Warnings

The IOL must be implanted in accordance to the following instructions for use. Improper use may impose risk to the patient's health.

Do not use if product sterility or quality is thought to be compromised due to damaged packaging.

Do not attempt to modify or alter the device or any of its components.

Do not use after expiration date.

Do not use without viscoelastic solution.

Do not reuse. Reuse may impose serious risk to the patient's health.

Do not resterilize as this can produce undesirable side effects.

Do not implant the lens in the anterior chamber or in the ciliary sulcus.

PRECAUTIONS:

- Store at room temperature. Avoid temperatures below 5°C and above 45°C. Do not freeze.
- The LuxGood™ IOL and the viscoelastic solution used should be stored for at least 30 min at minimum 20°C before surgery.
- A high level of surgical skill is required to implant intraocular lenses. Surgeon should be fully trained before attempting to implant intraocular lenses.
- The IOL model, its power and expiration date should be verified before opening the protective packaging and before opening the individual sterile pouch/blister.
- Sterility of the device is guaranteed only if the individual sterile pouch/blister has not been opened or damaged.
- Handle the lens carefully to avoid surface or haptics damages.
- The IOL must be carefully rinsed prior to implantation. Do not soak or rinse the IOL with any solution other than sterile balanced salt solution.
- Care should be taken to remove all viscoelastic from the anterior and posterior surfaces of the lens after implantation in the capsular bag.
- Care should be taken to verify that the lens is fully implanted into the capsular bag.
- If a YAG laser posterior capsulotomy is performed, assure that the laser beam is focused slightly behind the posterior capsule.

Ophthalmic kit for cataract surgery LUXGOOD™ YELLOW intraocular lens	Instructions for use for the Russian Federation version 2
Cutting Edge SAS, France	770, Rue Alfred Nobel, Immeuble le Nobel, 34000 Montpellier, France

PACKAGING / STERILIZATION:

The LuxGood™ IOL is individually package preloaded in an Accuject Pro delivery system in a blister. All models are supplied dry and should be opened under sterile conditions. A patient card and self-adhesive labels are supplied in the carton to provide traceability of the IOL. All models are sterilized by Ethylene Oxide.

DIRECTIONS FOR USE:

1. The physician should determine preoperatively the power of the lens to be implanted, using calculation methods. The recommended A-constants listed on the lens carton are guidelines only.
2. Prior to implanting, examine the lens package for model, power, expiration date and proper configuration.
3. Open the carton and inspect the sterilized blister containing the lens and the injection system. Make sure it is not damaged. Verify that the information on the blister label is consistent with the information on the outer box and the self-adhesive labels. Peel open the blister and transfer the device onto the sterile field.

4. Inspect the device and make sure it is intact. The lens should be rinsed thoroughly by flushing balanced salt solution through the tip of the cartridge (**FIG. 1**) and from the end of the loading chamber (**FIG. 2**). Wait for 30 seconds.

Please refer to the enclosed figures to fully prepare the device before implantation, taking care not to touch the IOL:

FIG. 3: Apply a generous amount of viscoelastic through the tip of the cartridge. Ensure that the viscoelastic fully fills the injector tip. Wait for 30 seconds.

FIG. 4: Close the wings of the loading chamber. When the "click-lock" mechanism is engaged, the lens is securely in place.

The lens is ready for implantation.

5. Inject the lens in the shortest possible time after fully preparing the lens.
6. Insert the cartridge tip in the incision with the bevel of the injector tip facing down to the patient's eye. Proceed by pushing the plunger forward smoothly but continuously until the lens is delivered into the capsular bag.
7. Discard the injection device.
8. The lens is to be implanted with the anterior side of the lens facing up towards the anterior side of the eye. The orientation of the IOL can be verified by visual inspection of the haptics. As

Ophthalmic kit for cataract surgery
LUXGOOD™ YELLOW intraocular lens

Instructions for use
for the Russian Federation
version 2

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illustrated in FIG. 5, when the haptic insertions are top right (A) and bottom left (B), you are looking at the anterior side of the lens.


FIG. 1

FIG. 2


FIG. 3

FIG. 4


FIG. 5


FIG. 6

Only viscoelastic solution of low to medium viscosity should be used. Only the injection system provided with the preloaded LuxGood™ lens should be used.

There are various surgical procedures that can be utilized, and the surgeon should select a procedure that is appropriate for the patient.

Any serious incident involving the device should be notified to the manufacturer and the competent authority of the Member State where the user is established.

RECYCLING:

Disposal of the injection system should be done in compliance with the local regulation of waste from health-care activities characterized by risk of infection.

After use, the injector for IOL implantation is classified as waste class B and must be disposed of in accordance with SanPin 2.1.3684-21.

Ophthalmic kit for cataract surgery LUXGOOD™ YELLOW intraocular lens	Instructions for use for the Russian Federation version 2
Cutting Edge SAS, France	770, Rue Alfred Nobel, Immeuble le Nobel, 34000 Montpellier, France

STORAGE AND TRANSPORTATION CONDITIONS:

Transportation conditions:

Temperature range from -10°C to +40°C and relative humidity from 85% to 90%.

Storage conditions:

Store at room temperature. Avoid temperatures below 5°C and above 45°C. Do not freeze.

EXPIRY DATE:

3 years.

LIFETIME:

IOL is intended to be permanently implanted in the body.

SYMBOLS:

YGOOD – Model name

D – Diopter

11 mm – Overall diameter

6 mm – Optic diameter

ϕ **T** – Overall diameter

ϕ **B** – Optic diameter

Serial number

Use before

Catalog number

Date of manufacture

Lot

Manufacturer

A-Constant (est)

Ophthalmic kit for cataract surgery LUXGOOD™ YELLOW intraocular lens	Instructions for use for the Russian Federation version 2
Cutting Edge SAS, France	770, Rue Alfred Nobel, Immeuble le Nobel, 34000 Montpellier, France

Ethylene oxide sterilization

Do not use if the package is damaged

Do not reuse

Do not resterilize

Storage temperature from +5°C

You should read the instructions for use before use

CE with the number of the inspection institution - confirms compliance with EU regulations

Carefully! Refer to instructions for use

Cutting Edge SAS, 770, Rue Alfred Nobel, Immeuble le Nobel, 34000 Montpellier, France

ORGANIZATION AUTHORIZED TO ACCEPT CLAIMS FROM CUSTOMERS:

Bausch Health LLC, Shabolovka St. 31, Bld. 5, Moscow, 115162, Russia,
Tel.: +7-495-510-28-79.

Specifications of Accuject Pro

Inner diameter at the head of the Cartridge tip	1.38 ± 0,05 mm
Outer tip diameter at the head of the Cartridge tip	1.66 ± 0.05 mm
Angle at the head of the tip	45 ± 2°
Total length	151 ± 0.1 mm
Length of the Loading Chamber	18,5 ± 1.0 mm
Width of the Loading chamber (1)	13,53 ± 0.04 mm
Length of the Handpiece	95.77 ± 0.05 mm

Ophthalmic kit for cataract surgery LUXGOOD™ YELLOW intraocular lens	Instructions for use for the Russian Federation version 2
Cutting Edge SAS, France	770, Rue Alfred Nobel, Immeuble le Nobel, 34000 Montpellier, France

Width of the Handpiece (1)	35.025 ± 0.03 mm
Width of the Handpiece (2)	8.046 ± 0.03 mm
Length of the Plunger	106.88 ± 0.05 mm
Width of the Plunger (1)	14.927 ± 0.07 mm
Width of the Plunger (2)	5.483 ± 0.03 mm
Width of the Plunger (3)	1.755 ± 0.03 mm
Length of the Silicone tip	10.05 ± 0.05 mm
Width of the Silicone tip (1)	2.348 ± 0.04 mm
Width of the Silicone tip (2)	3.532 ± 0.04 mm
Length of the Cartridge tip	29.11 ± 0.08 mm
Width of the Cartridge tip (1)	10.231 ± 0.04 mm
Width of the Cartridge tip (2)	5.537 ± 0.04 mm
Pressing force on the plunger Injection force	20N ± 2,2 N

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